

SPECIFICATION & DATA SHEET

AREA VALVE SERVICE UNIT (AVSU)

The MIPS Area Valve Service Unit (AVSU) is designed for economy, maximum design life and reliability whilst requiring minimum maintenance. They are available as a standard product with the following design features:

- CE marked
- Powder-coated steel enclosure
- Fully removable quick-release door for easy access
- Reversible door hinge (left or right)
- Complies fully with the requirements of HTM 02-01
- All valves fitted with BS 5682 NIST connectors
- Sealed shut-off NIST valves
- All valves fitted with flat-faced 'o' ring unions
- Standard pipework sizes 15-28mm
- Each valve supplied with blanking facility
- Every AVSU fully tested prior to despatch
- Every valve tested for tightness (no pressure pass when closed)



In addition to all of the above, MIPS can offer the following options:

- Flush surround and disposable plaster shield
- 35-54mm pipework valves
- Back or top entry
- Suited key locks