

SPECIFICATION & DATA SHEET

MEDICAL AIR PLANT

The MIPS Air Plant provides a centralised source of medical air. The systems are designed for economy, reliability and durability, whilst offering the following features as standard;

- Fully complies with HTM 02-01
- Compact and aesthetically designed
- All inter-connecting pipework and wiring is included
- Flexible connections on compressor pipe-work
- Anti-vibration mounts on compressor units
- Fail-safe control system in event of power loss
- Each machine fitted with adjacent electrical isolator
- Control panel fitted with highly visible indicators
- Integral alarm system interface indicator
- Control includes full operating and indicating system
- Alarm and BMS outputs for remote status indication
- Line fault monitoring of standby pressure sensor wiring
- 400kpa and 700kpa output versions available
- All pressure switches and regulators are factory set for ease of installation
- All plants are tested and certified before despatch
- Each plant is covered by a full warranty



Duplex Medical Air Plant



Triplex Medical Air Plant

In addition to all of the above, MIPS can custom build air plants to your specific requirements. Here are some examples available;

- Acoustic hooded machines
- Special layouts
- Screw compressors
- Other supply voltages
- Combination breathing/tool air plant