

SPECIFICATION & DATA SHEET

MEDICAL AUTOMATIC MANIFOLD

The MIPS Medical Automatic Manifold is designed for economy, reliability and durability, while offering the following features as standard;

- Fully complies with HTM 02-01
- Fully automatic
- Two-stage regulators for good flow characteristics
- Robust wall-mounted framework
- Copper pipeline union connections flat-faced with 'o' ring
- Pre-set regulators and pressure sensors
- Highly visible gauges for cylinder contents and pipeline pressure
- Tailpipe check-valves and integral non-return valves
- Cylinder support chains and racks included
- Includes safety relief valve and isolation valve
- Compact and aesthetically designed enclosure
- Removable cover to allow easy maintenance
- Fail-safe control system in the event of a power loss
- Control panel fitted with highly visible indicators
- Integral alarm system interface indicator
- Alarm outputs for remote status indication
- Line fault monitoring of standby pressure sensor wiring
- 400kpa, 700kpa and 1050kpa output versions available
- All manifolds are tested and certified before despatch
- Each manifold is covered by a full warranty



The right is reserved to change the specification of the items described herein at any time.