

SPECIFICATION & DATA SHEET

MEDICAL GAS ALARM: SDX15

The MIPS Medical Gas Alarm provides a centralised monitor for up to five gas systems. The systems are designed for economy, reliability and durability, whilst offering the following features as standard;

- Fully complies with HTM 02-01
- Up to 5 gases can be included
- Compact and aesthetically designed
- 24-hour battery back-up included in the event of power loss
- Integral alarm system interface indicator
- BMS outputs for remote status indication
- All alarms are fully tested and certified before despatch
- Each alarm is covered by a full warranty

