

SPECIFICATION & DATA SHEET

MEDICAL GAS TERMINAL UNIT

The MIPS Medical Gas Terminal Unit provides a secure anti-confusion supply of specific gas, at the correct pipeline pressure at point-of-use.

- Meets the requirements of HTM 02-01 and IS EN 737
- All components are indexed to ensure gas specificity
- Forms part of an integrated range including AGSS to BS 6483 and IS EN 737
- Hassle-free installation
- Manufactured to standards meeting the requirements of ISO 9001
- Efficient maintenance, removal of 2nd fix without shutting down the system
- Complete with copper pie stub complying to WKO 82/1
- Available in all formats; flush, surface, bedhead trunking, flexible pendant, multi-movement pendant

