


SPECIFICATION & DATA SHEET

OXYGEN FLOWMETERS FOR PIPELINE USAGE

MIPS oxygen flowmeters offer you a user-friendly, robust and accurate medical device whose calibration cannot be altered.

Maintain suitability for direct connection to the terminal unit using B.S. Probe. Remote connection can be achieved using an anti-static hose and B.S. Probe. They are available as a standard product with the following design features:

- CE marked
- Highly shock-resistant, polycarbonate reading scale and scale cover
- Flow rate: 15 l/min
- Filter positioned at the scale inlet, ensuring patient and flowmeter protection
- Extremely robust chromium-plated brass body
- Easy to maintain knob offers optimum adjustment quality
- Vertical positioning
- Tubing nipple
- Every unit rigorously tested prior to despatch
- Individual serial number stamped on the body of each flowmeter for identification purposes
- Complies fully with the requirements of HTM 02-01



In addition to all of the above, MIPS can offer the following options:

- Twin flowmeters
- Humidifiers & adaptors
- Oxygen tube & mask
- Spare tubing nipples

The right is reserved to change the specification of the items described herein at any time.