

SPECIFICATION & DATA SHEET

PORTABLE OR MOBILE SUCTION TROLLEY

The MIPS Suction Trolley is designed for economy, maximum design life and reliability, whilst requiring minimum maintenance.

- CE marked
- Powder-coated steel frame
- Totally transparent, ensuring perfect visibility of the liquid contents
- Silicone ring pressure-fit lid available
- Environmentally-friendly, PCV-free suction liners
- Single or twin versions of RBTM3 vacuum regulator or twin RVTM3 available
- Complies fully with requirements of HTM 02-01
- Every unit rigorously tested prior to despatch
- Batch number marked on the base of each jar for identification
- 3m long silicone hose included
- 2L collection jars available in polysulfone
- 3m long connection to the central vacuum source with remote probe
- One twin catheter holder



The right is reserved to change the specification of the items described herein at any time.