

SPECIFICATION & DATA SHEET

RIGID PENDANTS

There are three basic types of MIPS Rigid Pendants:

Mini – available with up to 4 gases and 4 single switched sockets

Midi – available with up to 9 gases and 8 double switched sockets

Maxi – available with up to 9 gases and 8 double switched sockets with adjustable length by 300mm

Each is designed for aesthetics, economy, maximum design life and reliability, whilst requiring minimum maintenance. They are available as standard products with the following design features:

- Comply fully with the requirements of HTM 02-01
- Sturdy monocoque design
- Manufactured from steel and epoxy powder coated
- Standard colour white
- All gases fitted with BS 5682 NIST hoses for easy installation
- First fix available separately and supplied with blanking facility
- Every rigid pendant is pressure tested prior to despatch
- Every rigid pendant is electrically tested prior to despatch



In addition to all of the above, MIPS can **custom build** rigid pendants to your specific requirements. Here are some examples:

- Other colours are available
- Rigid copper pipework
- Special electrical connections e.g. data, video etc
- Other makes of gas outlets are available

The right is reserved to change the specification of the items described herein at any time.