

## SPECIFICATION & DATA SHEET

### SUCTION CONTROLLER & COLLECTION JAR

The MIPS Suction Controller and Collection Jar is designed for economy, maximum design life, stability and reliability, whilst requiring minimum maintenance. The suction regulator is suitable for direct connection to a vacuum terminal unit using DIN probe. You can also complete a remote connection using additional vacuum flexible hose and DIN probe.

They are available as a standard product with the following design features:

- CE Marked
- Anti-overflow safety device ensuring ease of maintenance
- Totally transparent, ensuring perfect visibility of the liquid contents
- 150cc polycarbonate safety jar with an antibacterial filter
- Silicone ring lid on 2 litre collection jar
- Easy to use
- Environmentally-friendly, PVC-free suction liners
- Fully complies with the requirements of HTM 02-01
- Every unit rigorously tested prior to despatch
- Batch number marked on the base of each jar for identification



In addition to all of the above, MIPS can offer the following options:

- Silicone suction hoses (25cm coil)
- Spare antibacterial filters
- Wide range of accessories (rail adaptors, wall brackets etc)
- Disposable suction liners